

Contents

BioPipelines User Manual	1
Table of Contents	1
Architecture Overview	1
Core Concepts	2
Tool Abstract	2
Structure Generation Tools	3
Sequence Design Tools	6
Compound Generation Tools	8
Analysis Tools	9
Filtering Tools	10
Optimization Tools	11
Utility Tools	12
Pipeline Examples	13
Troubleshooting	17

BioPipelines User Manual

Table of Contents

1. Architecture Overview
2. Core Concepts
3. Tool Abstract
4. Structure Generation Tools
5. Sequence Design Tools
6. Compound Generation Tools
7. Analysis Tools
8. Filtering Tools
9. Optimization Tools
10. Utility Tools
11. Pipeline Examples
12. Troubleshooting

Architecture Overview

BioPipelines is a system that **generates bash scripts** for bioinformatics workflows. The pipeline itself does not execute computations directly - instead, it predicts and prepares the filesystem structure that will exist at SLURM runtime, then generates bash scripts to be executed when the SLURM job runs.

Execution Model

- **Pipeline Role:** Generates bash scripts and predicts filesystem structure
- **Execution Time:** Scripts are executed only at SLURM runtime
- **HelpScripts:** Utility scripts designed to be executed during SLURM job execution

- **Filesystem Prediction:** Pipeline pre-determines output paths and directory structure before job submission

Key Directories

- **PipelineScripts/:** Tool classes that generate pipeline configurations
 - **HelpScripts/:** Runtime utility scripts executed during SLURM jobs
 - **Output Structure:** Each tool creates numbered directories (1_RFdiffusion/, 2_ProteinMPNN/, etc.)
-

Core Concepts

Datasheets

All tools output standardized CSV datasheets that serve as the primary interface between pipeline steps:

- **Structure Generation:** id, structure_file, [tool-specific columns]
- **Sequence Design:** id, sequence, source_id, [metrics]
- **Compound Generation:** id, format, smiles, ccd
- **Analysis Results:** id, source_structure, [metric_columns]

Tool Chaining

Tools are connected through their standardized outputs:

```
rfd = RFdiffusion(...)
pmpnn = ProteinMPNN(input=rfd.output)
alphafold = AlphaFold(input=pmpnn.output)
```

Environment Management

Each tool specifies its conda environment requirements for SLURM execution.

Tool Abstract

Tool Structure

Every BioPipelines tool follows this abstract pattern:

```
class ToolName(BaseConfig):
    TOOL_NAME = "ToolName"

    def __init__(self, input=None, **parameters):
        # Tool-specific initialization

    def get_output_files(self) -> Dict[str, List[str]]:
        # Returns predicted output file paths
        return {
```

```

        "structures": [...],      # PDB/CIF files
        "sequences": [...],      # FASTA files
        "compounds": [...],      # Compound files
        "datasheets": {...},     # CSV datasheets
        "output_folder": self.output_folder
    }

    @property
    def output(self) -> StandardizedOutput:
        # Standardized output interface

```

Input/Output Pattern

- **Input:** input=previous_tool.output or direct parameters
- **Output:** Standardized datasheets + files accessible via tool.output
- **Chaining:** Tools consume the output property of upstream tools

Datasheet Standardization

- **ID Column:** Always present, unique identifier for each item
- **Source Tracking:** source_id, source_structure columns link to upstream tools
- **File References:** Absolute paths to generated files
- **Metrics:** Tool-specific analysis columns

Structure Generation Tools

RFdiffusion

Purpose: Generate novel protein structures using diffusion models

Usage

```

rfd = RFdiffusion(
    name="my_designs",
    pdb="path/to/template.pdb",      # Optional template
    contigs="A1-100/50-80",          # Chain definitions
    length="200-300",                # Target length range
    num_designs=10,                  # Number of structures
    environment="rfdiffusion"
)

```

Parameters

- **name:** Job identifier
- **pdb:** Template PDB file (optional)
- **contigs:** Chain length specifications
- **length:** Total length or range
- **num_designs:** Number of structures to generate
- **environment:** Conda environment name

Input

- Optional template PDB file
- Contig and length specifications

Output Datasheet: `rfdiffusion_results.csv` | Column | Description | |-----|-----| |
| `id` | Structure identifier | | `structure_file` | Path to generated PDB | | `fixed` | Fixed regions
(PyMOL selection) | | `designed` | Designed regions (PyMOL selection) | | `exists` | Boolean, file
existence check |

Files: PDB structures in output directory

AlphaFold

Purpose: Predict protein structures from sequences using AlphaFold

Usage

```
af = AlphaFold(  
    input=pmpnn.output,           # From ProteinMPNN  
    name="predictions",  
    rank=True,                    # Rank by confidence  
    num_relax=3,                  # Number to relax  
    environment="alphafold"  
)
```

Parameters

- `input`: Sequences from upstream tool
- `name`: Job identifier
- `num_relax`: Number of best models to relax
- `num_recycle`: Recycling iterations (default 3)
- `environment`: Conda environment

Input

- Sequence datasheet from ProteinMPNN or similar
- FASTA files with protein sequences

Output Datasheet: Inherits from input sequences | Column | Description | |-----|-----| |
| `id` | Sequence identifier | | `source_id` | Original sequence ID | | `sequence` | Protein sequence | |
`structure_file` | Path to predicted PDB |

Files: PDB structures, confidence scores, MSAs

Boltz2

Purpose: Predict protein-ligand complex structures

Usage

```
boltz = Boltz2(  
    proteins=pmpnn.output,          # Protein sequences  
    ligands=compound_lib.output,    # Compound library  
    name="complexes",  
    num_samples=5,                  # Sampling iterations  
    environment="boltz2"  
)
```

Parameters

- **proteins**: Protein sequences (FASTA or upstream tool)
- **ligands**: Compound library or SMILES
- **config**: Configuration files
- **msas**: Multiple sequence alignments
- **name**: Job identifier
- **num_samples**: Number of samples per prediction
- **environment**: Conda environment

Input

- Protein sequences (FASTA files)
- Ligand library (compound datasheets)
- Optional MSAs and configuration files

Output Multiple Datasheets:

1. **Confidence**: `boltz_confidence_scores.csv` | Column | Description | |———|—————| | **id** | Complex identifier | | **input_file** | Input configuration name | | **confidence_score** | Overall confidence | | **ptm** | Predicted Template Modeling score | | **iptm** | Interface PTM score |
2. **Sequences**: `boltz_sequences.csv` | Column | Description | |———|—————| | **id** | Sequence identifier | | **sequence** | Protein sequence |
3. **Compounds**: `boltz_compounds.csv` | Column | Description | |———|—————| | **id** | Compound identifier | | **smiles** | SMILES string | | **format** | Format type (SMILES/CCD) | | **ccd** | CCD identifier if applicable |

Files: CIF structure files, JSON confidence scores

RFdiffusionAllAtom

Purpose: Generate ligand-aware protein structures with all-atom diffusion

Usage

```
rfd_aa = RFdiffusionAllAtom(  
    ligand="ZIT",                  # Ligand identifier
```

```

pdb="template.pdb",           # Template structure
contigs="A1-100,10-20",      # Contig specification
num_designs=5,               # Number of designs
active_site=True,            # Use active site model
environment="rfdiffusion"
)

```

Parameters

- **ligand**: Ligand identifier (e.g., ‘ZIT’, ‘RFP’)
- **pdb**: Input PDB template (optional)
- **contigs**: Contig specification for design regions
- **inpaint**: Inpainting regions specification
- **num_designs**: Number of structures to generate
- **active_site**: Use active site model for small motifs
- **steps**: Diffusion steps (default 200)
- **ppi_design**: Enable protein-protein interaction design
- **ppi_hotspot_residues**: List of hotspot residues for PPI
- **environment**: Conda environment

Input

- Template PDB file with ligand context
- Ligand specifications and design regions
- Optional hotspot residue definitions

Output Datasheet: `rfdiffusion_allatom_results.csv` | Column | Description |

id	Structure identifier	source_id	Template source ID	structure_file	Path to generated PDB	fixed	Fixed regions (PyMOL selection)	designed	Designed regions (PyMOL selection)	contigs	Contig specification used	time	Generation time	status	Completion status
----	----------------------	-----------	--------------------	----------------	-----------------------	-------	---------------------------------	----------	------------------------------------	---------	---------------------------	------	-----------------	--------	-------------------

Files: PDB structures with ligand contexts

Sequence Design Tools

ProteinMPNN

Purpose: Design protein sequences for given structures

Usage

```

pmpnn = ProteinMPNN(
    input=rfd.output,           # Structure input
    name="sequences",
    num_sequences=50,           # Sequences per structure
    temperature=0.1,           # Sampling temperature
    environment="proteinmpnn"
)

```

Parameters

- **input**: Structures from upstream tool
- **name**: Job identifier
- **num_sequences**: Number of sequences per structure
- **temperature**: Sampling temperature
- **batch_size**: Processing batch size
- **environment**: Conda environment

Input

- Structure datasheet from RFdiffusion or similar
- PDB files with protein structures

Output Datasheet: `pmpnn_results.csv` | Column | Description | |
Sequence identifier | | `source_id` | Source structure ID | | `source_structure` | Path to input PDB |
| `sequence` | Designed protein sequence | | `score` | ProteinMPNN score | | `seq_recovery` | Sequence recovery rate |

Files: FASTA files with designed sequences

LigandMPNN

Purpose: Design protein sequences considering ligand interactions

Usage

```
lmpnn = LigandMPNN(  
    input=rfd.output,           # Structure input  
    ligand_params="ligand.json", # Ligand parameters  
    name="ligand_designs",  
    num_sequences=30,  
    environment="ligandmpnn"  
)
```

Parameters

- **input**: Structures with ligands
- **ligand_params**: Ligand parameter file
- **name**: Job identifier
- **num_sequences**: Sequences per structure
- **temperature**: Sampling temperature
- **environment**: Conda environment
- **fixed**: Selection of positions to keep fixed
- **designed**: Selection of positions to redesign

Input

- Protein structures with bound ligands

- Ligand parameter configuration

Output Datasheet: `lmpnn_results.csv` | Column | Description | |———|—————| | **id** | Sequence identifier | | **source_id** | Source structure ID | | **sequence** | Designed sequence | | **ligand_binding_score** | Binding affinity score |

Files: FASTA files with ligand-aware sequences

Compound Generation Tools

CompoundLibrary

Purpose: Generate and process compound libraries for structure prediction

Usage

```
compounds = CompoundLibrary(
    library={"core": "C1CCCCC1", "R1": ["F", "Cl", "Br"]},
    name="my_library",
    primary_key="core",
    max_compounds=1000,
    environment="rdkit"
)
```

Parameters

- **library:** Dictionary or file paths with compound definitions
- **name:** Library identifier
- **primary_key:** Primary library component
- **max_compounds:** Maximum number to generate
- **covalent:** Generate covalent variants
- **environment:** Conda environment

Input

- **Dictionary:** {"scaffold": "SMILES", "R1": ["substitutions"]}
- **File Paths:** CSV/TXT files with SMILES
- **SMILES String:** Single compound

Output Datasheet: `{name}_compounds.csv` | Column | Description | |———|—————| | **id** | Compound identifier | | **smiles** | SMILES string | | **format** | Format type (SMILES) | | **ccd** | CCD code (if applicable) |

Optional Datasheet: `{name}_compound_properties.csv` | Column | Description | |———|—————| | **id** | Compound identifier | | **MW** | Molecular weight | | **LogP** | Lipophilicity | | **HBD** | Hydrogen bond donors | | **HBA** | Hydrogen bond acceptors |

Files: JSON library definitions, property calculations

Analysis Tools

ResidueAtomDistance

Purpose: Calculate distances between residue selections and atom coordinates

Usage

```
distance = ResidueAtomDistance(  
    input=boltz.output,          # Structures to analyze  
    residue="D in IGDWK",  
    atom_coords="LIG.C1",  
    name="active_site_distance" # Metric name  
)
```

Parameters

- **input:** Structures from upstream tool
- **residue_selection:** PyMOL-style residue selection
- **atom_coords:** Target atom coordinates [x, y, z]
- **name:** Analysis identifier
- **distance_cutoff:** Maximum distance threshold

Input

- Structure datasheet with PDB files
- Selection criteria and target coordinates

Output Datasheet: {name}_distances.csv | Column | Description | | id |
Structure identifier | | source_structure | Path to analyzed PDB | | {metric_name} | Distance
measurement |

Confidence

Purpose: Extract confidence scores (pLDDT) from protein structures

Usage

```
confidence = Confidence(  
    input=alphafold.output,      # Structures to analyze  
    name="plddt"                 # Custom metric name  
)
```

Parameters

- **input:** Structure datasheet from prediction tools
- **name:** Analysis identifier
- **metric_name:** Custom name for confidence metric column
- **residue_range:** Specific residues to analyze (optional)

Input

- Structure datasheet with PDB files
- Protein structures with B-factor confidence scores

Output Datasheet: {name}_confidence.csv | Column | Description | |———|—————| | id |
Structure identifier | | source_structure | Path to analyzed PDB | | <name> | pLDDT value |

Confidence Score Interpretation

- **Very high confidence:** >90 (highly accurate)
 - **High confidence:** 70-90 (generally accurate)
 - **Low confidence:** 50-70 (potentially inaccurate)
 - **Very low confidence:** <50 (likely inaccurate)
-

Filtering Tools

Filter

Purpose: Filter datasheets using pandas query expressions

Usage

```
filter_tool = Filter(  
    input=analysis.output,          # Datasheet to filter  
    expression="plddt_avg > 70 & contact_count > 50",  
    name="high_quality",  
    max_items=100                  # Limit results  
)
```

Parameters

- **input:** Datasheet from analysis tools
- **expression:** Pandas query expression
- **name:** Filter identifier
- **combination:** “AND” or “OR” for multiple criteria
- **max_items:** Maximum items to keep
- **score_weights:** Weighting for scoring criteria

Input

- Analysis datasheet with metrics
- Filter expression criteria

Output Datasheet: Filtered version of input - Same columns as input - Only rows meeting filter criteria - Additional metadata about filtering

Expression Examples

```
# Simple filters
"confidence > 0.8"
"distance < 5.0"

# Complex filters
"(plddt_avg > 70) & (contact_count > 50)"
"score > 0.9 | (confidence > 0.8 & distance < 3.0)"

# String matching
"id.str.contains('design')"
```

Optimization Tools

SelectBest

Purpose: Select the single best item from analysis results for iterative optimization

Usage

```
best = SelectBest(
    input=analysis.output,      # Analysis results
    metric="binding_affinity",  # Primary optimization metric
    mode="max",                 # Maximize or minimize
    name="best_design"
)
```

Parameters

- **input:** Analysis datasheet with metrics
- **metric:** Primary metric to optimize
- **mode:** “max” or “min” to maximize or minimize metric
- **weights:** Dict of {metric_name: weight} for multi-objective selection
- **tie_breaker:** How to handle ties (“first”, “random”, or metric name)
- **composite_function:** How to combine metrics (“weighted_sum”, “product”, “min”, “max”)
- **name:** Name for the selected best item

Input

- Analysis datasheet with quantitative metrics
- Selection criteria and optimization parameters

Output Datasheet: Single row with best selected item - Same columns as input datasheet - Only the single best item based on criteria - Additional metadata about selection process

Multi-Objective Selection

```
# Weighted combination of multiple metrics
best_multi = SelectBest(
    input=combined_analysis,
    metric="composite_score",
    weights={"binding_affinity": 0.6, "plddt_avg": 0.4},
    mode="max",
    tie_breaker="binding_affinity"
)
```

Use Cases

- **Iterative Design:** Select best design for next optimization cycle
- **Multi-Objective:** Balance multiple competing objectives
- **Quality Control:** Pick highest quality result from batch
- **Resource Optimization:** Select most promising candidate for expensive analysis

Utility Tools

MergeDatasheets

Purpose: Combine multiple datasheets with intelligent column handling

Usage

```
merger = MergeDatasheets(
    input_datasheets=[analysis1.output, analysis2.output],
    join_key="id",
    name="combined_analysis"
)
```

Parameters

- **input_datasheets:** List of datasheet sources
- **join_key:** Column to merge on (default: "id")
- **name:** Merger identifier
- **handle_duplicates:** Strategy for duplicate columns

Input

- Multiple datasheets with compatible ID columns
- Join specifications

Output Datasheet: Combined datasheet - All unique columns from input datasheets - Outer join preserving all data - Collision handling for duplicate column names

LoadOutput

Purpose: Load results from completed pipeline runs

Usage

```
loaded = LoadOutput(  
    result_file="output.json",    # Pipeline result metadata  
    validate_files=True  
)
```

Parameters

- `result_file`: Pipeline output JSON file
- `validate_files`: Check file existence
- `load_datasheets`: Load datasheet contents

Input

- Pipeline result JSON metadata
- Output directory structure

Output

- Access to all pipeline outputs
 - Datasheet loading and validation
 - File existence verification
-

Pipeline Examples

Basic Protein Design Pipeline

```
from PipelineScripts import *  
  
# Generate novel structures  
rfd = RFdiffusion(  
    name="novel_designs",  
    contigs=["A1-150"],  
    num_designs=20,  
    environment="rfdiffusion"  
)  
  
# Design sequences  
pmpnn = ProteinMPNN(  
    input=rfd.output,  
    name="sequences",  
    num_sequences=10,  
    temperature=0.1,  
    environment="proteinmpnn"
```

```

)

# Predict structures
alphafold = AlphaFold(
    input=pmpnn.output,
    name="predictions",
    rank=True,
    num_relax=5,
    environment="alphafold"
)

# Analyze quality
analysis = Confidence(
    input=alphafold.output,
    name="quality"
)

# Filter best results
filter_best = Filter(
    input=analysis.output,
    expression="plddt_avg > 80 & contact_count > 100",
    name="high_quality",
    max_items=5
)

# Create pipeline
pipeline = Pipeline([rfd, pmpnn, alphafold, analysis, filter_best])
pipeline.save()

```

Protein-Ligand Complex Design

```

# Generate compound library
compounds = CompoundLibrary(
    library={"scaffold": "c1ccccc1", "R1": ["F", "Cl", "N"]},
    name="drug_library",
    max_compounds=50,
    environment="rdkit"
)

# Design protein structures
rfd = RFdiffusion(
    name="binding_sites",
    length="100-200",
    num_designs=10,
    environment="rfdiffusion"
)

# Predict complexes

```

```

boltz = Boltz2(
    proteins=rfd.output,
    ligands=compounds.output,
    name="complexes",
    num_samples=3,
    environment="boltz2"
)

# Analyze binding
distance_analysis = ResidueAtomDistance(
    input=boltz.output,
    residue_selection="resi 50-80",
    atom_coords=[0.0, 0.0, 0.0],
    name="binding_distance"
)

# Filter viable complexes
viable = Filter(
    input=distance_analysis.output,
    expression="binding_distance < 4.0 & confidence_score > 0.7",
    name="viable_complexes"
)

pipeline = Pipeline([compounds, rfd, boltz, distance_analysis, viable])

```

Iterative Design Optimization

```

# Initial design
rfd = RFdiffusion(
    name="initial_designs",
    contigs=["A1-100"],
    num_designs=50,
    environment="rfdiffusion"
)

# Analyze quality
confidence = Confidence(
    input=rfd.output,
    name="initial_confidence"
)

# Select best for refinement
best_initial = SelectBest(
    input=confidence.output,
    metric="plddt_avg",
    mode="max",
    name="best_initial"
)

```

```

# Design sequences for best structure
pmpnn = ProteinMPNN(
    input=best_initial.output,
    name="refined_sequences",
    num_sequences=20,
    environment="proteinmpnn"
)

# Predict refined structures
alphafold = AlphaFold(
    input=pmpnn.output,
    name="refined_structures",
    environment="alphafold"
)

# Multi-objective analysis
final_confidence = Confidence(
    input=alphafold.output,
    name="final_confidence"
)

binding_analysis = ResidueAtomDistance(
    input=alphafold.output,
    residue_selection="resi 50-70",
    atom_coords=[0.0, 0.0, 0.0],
    name="binding_site"
)

# Combine metrics
combined = MergeDatasheets(
    input_datasheets=[final_confidence.output, binding_analysis.output],
    name="combined_metrics"
)

# Select optimal design
final_best = SelectBest(
    input=combined.output,
    metric="composite_score",
    weights={"plddt_avg": 0.7, "binding_site_distance": 0.3},
    mode="max",
    name="optimal_design"
)

pipeline = Pipeline([rfd, confidence, best_initial, pmpnn, alphafold,
                    final_confidence, binding_analysis, combined, final_best])

```


Multi-Analysis Filtering Pipeline

```
# Load existing structures
structures = LoadOutput(result_file="previous_run.json")

# Multiple analyses
confidence_analysis = Confidence(
    input=structures.output,
    name="confidence"
)

distance_analysis = ResidueAtomDistance(
    input=structures.output,
    residue_selection="resi 100-120",
    atom_coords=[10.0, 15.0, 20.0],
    name="active_site"
)

# Combine analyses
combined = MergeDatasheets(
    input_datasheets=[confidence_analysis.output, distance_analysis.output],
    name="all_metrics"
)

# Complex filtering
final_filter = Filter(
    input=combined.output,
    expression="(plddt_avg > 75) & (active_site_distance < 5.0) & (helix_content > 0.3)",
    name="final_selection",
    max_items=10
)
```

Troubleshooting

Common Issues

- 1. File Path Errors** **Problem:** Tools cannot find input files **Solution:** Ensure upstream tools completed successfully and files exist
- 2. Environment Conflicts** **Problem:** Tool fails due to conda environment issues **Solution:** Verify environment names match available conda environments
- 3. Memory Issues with Large Libraries** **Problem:** Compound libraries too large for processing **Solution:** Use `max_compounds` parameter to limit library size
- 4. Filter Expression Errors** **Problem:** Filter expressions fail with syntax errors **Solution:** Use pandas query syntax, quote string columns

5. Missing Dependencies **Problem:** Tools fail due to missing software dependencies **Solution:** Check that all required tools are installed in specified environments

Debugging Tips

1. **Check Log Files:** Each tool generates logs in the main pipeline directory
2. **Validate Datasheets:** Inspect CSV outputs to understand data flow
3. **Use Small Test Cases:** Start with minimal examples before scaling up
4. **File Existence:** Verify all input files exist and are accessible
5. **Environment Testing:** Test tools individually in their target environments

Best Practices

1. **Incremental Development:** Build pipelines step by step
2. **Data Validation:** Check intermediate outputs before chaining tools
3. **Resource Planning:** Consider computational requirements for each tool
4. **Error Handling:** Plan for failed jobs and restart strategies
5. **Documentation:** Document custom parameters and pipeline logic